

HOW TO ACCESS A CUSTOMIZABLE REDCAP EXAMPLE FOR A TRACKING AND MONITORING SYSTEM

1. First, you have access to REDCap. If yes, go to the REDCap Shared Library.
2. Access the REDcap Shared Library option either on the Project Setup Tab -> "Design Data Collection instruments..." section, or via the "Import" option in the Online Designer.
3. Search for a keyword from the form in the library. The 3 instruments from the Transition Registry are Patient Demographics, Pediatric Transition Registry and Adult Transition Registry.
4. Open the dropdown for the located document, and select "Import into My REDCap project."